1. GENERAL INFORMATION

Intended Use: IV Administration Set is intended for the infusion, transfer, and/or administration of medications, drugs or fluids.

Indications: Patient medical condition in needs of infusion.

Contraindications: Not known

Medical personnel responsibility for patient safety: In order to competently manage a patient the medical personnel is required to:

- respond appropriately (adjust treatments as necessary)
- observe for complications and troubleshoot as needed

Warning notices and symbols: The following symbols are used for important information in the Instruction for use, on the packaging and on the product:

	Caution, consult accompanying documents		Use by date
	Manufactured by		Sterilized using ethylene oxide
	CE mark of compliance		Do not use if package is damaged
	Handling mark on package – FRAGILE		For single use only
	Handling mark on package – KEEP DRY		Recyclable material
	Handling mark on package – THIS SIDE UP		Consult instruction for use
	Do not re-sterilise		Non latex
	PHT free		

2. DEVICE DESCRIPTION

IV Administration Set may consist of various components: *Slide clamps, Needle free valve, Tubing, Luer Locks, Caps for Luer Locks.*

3. OPERATION

 Its highly recommended to consult instruction prior to use

1. Visually check packaging integrity
-  Product is sterile, do not use if package damaged or open
2. Using antiseptic technique, open the set from the sterile packing.
3. Ensure that all connections are secure, *as shown on Fig. 1.*

 All connections should be finger tightened. Over tightening can cause cracks and leaks to occur that could result in embolism and or exposure to biohazards.

4. Close the slide clamps on the IV set.
5. Remove protective covering from spike tip and take special attention that exposed end of drip chamber spike will not come into contact with anything, *per Fig. 2.*

 Exposed end of drip chamber spike should not come into contact with anything, otherwise it can be contaminated.

6. Connect saline bag to the drip chamber set.
7. Prime the tubing and remove air from the IV flush solution bag per hospital policy

 If all air is not removed from the tubing, air may be forced into the patient's vascular system and cause embolism.

8. Fill the drip chamber halfway with flush solution by squeezing the drip chamber, *as shown on Fig. 3.* Open the clamp.

Fig. 1

Fig. 2

Fig. 3

Check for fluid leakage before and during the procedure. Leaks can result in the loss of sterility, fluid and/or air embolism. If a product leaks before or during use, retighten the leaking connection or replace the product

9. Connect IV tubing to the catheter, as shown on Fig. 4. Flush system per hospital policy

Fig. 4

Operation with infusion syringe pump

1. Connect the IV administration set to the syringe.
2. Place the syringe in the infusion syringe pump, as shown on Fig. 5.
3. Use the infusion syringe pump in “fluid push” option.
4. Pressurize the IV flush solution bag. Flow rate will be managed by infusion syringe pump.
5. Connect IV tubing to the catheter.
6. Flush system per hospital policy.


Fig.5

4. STORAGE AND DISPOSAL:

Keep dry

Use by date

Do not re-sterilize, *Improper reuse of the medical device can lead to potentially spread life-threatening infections*

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