Instructions for use – Contrast Media Delivery Set

1. GENERAL INFORMATION

Intended Use:

Contrast Media Delivery Set is intended for the infusion and/or administration of contrast media.

Indications:

Is indicated for a patients that undergo computed tomography (contrast CT), projectional radiography, and fluoroscopy procedures , in aim to present internal organs by minimal invasive approach

Contraindications:

The associated risks can be attributed to cardiac and non-cardiac complications. Specific disease states pertaining to the general medical profile of the patient (older age, renal insufficiency, uncontrolled diabetes mellitus, and morbid obesity) can increase the risk for complications. The underlying vascular status of the patient can further predispose to adverse events.

Medical personnel responsibility for patient safety:

In order to competently manage a patient the medical personnel is required to:

- respond appropriately (adjust treatments as necessary)
- observe for complications and troubleshoot as needed

Warning notices and symbols:

The following symbols are used for important information in the Instruction for use, on the packaging and on the product:

	Caution, consult accompanying documents		Use by date
	Manufactured by		Sterilised using ethylene oxide
	CE mark of compliance		Do not use if package is damaged
	Handling mark on package – FRAGILE		For single use only
	Handling mark on package – KEEP DRY		Recyclable material
	Handling mark on package – THIS SIDE UP		Consult instruction for use
	Do not re-sterilise		Non latex
	PHT free		

2. DEVICE DESCRIPTION

Contrast Media delivery Set may consist of various components: *Tubing line, Roller clamp, Syringe, Luer Locks, Spike Set, High Pressure Manifold, Contrast Media Set, Dual Check Valve, Drip Chamber Set*

3. OPERATION

Its highly recommended to consult instruction prior to use.


1. Visually check packaging integrity.

Product is sterile, do not use if package damaged or open.

2. Using antiseptic technique, open the set from the sterile packing.
3. Ensure that all connections are secure, *as shown on Fig. 1.*

Fig. 1

All connections should be finger tightened. Over tightening can cause cracks and leaks to occur that could result in embolism and or exposure to biohazards.

4. Remove protective covering from spike tip and take special attention that exposed end of drip chamber spike will not come into contact with anything. Repeat same procedure with Drip Chamber set spike. (see Figure 2)


Fig. 2

Exposed end of drip chamber spike should not come into contact with anything, otherwise it can be contaminated.

5. Connect contrast line to spike set, saline bag to the fluid reservoir drip chamber set (see Figure 3).

6. Connect syringe to dual check valve through Luer lock connection.
7. Clear the contrast line and tubing line with saline by turning the stopcock on manifold unit so the both lines are opened to syringe and locked to the patient line.

Fig. 3

8. Evacuate the air from the contrast media set and tubing line with saline by using of syringe.
9. Evacuate air from the syringe by closing the stop cock on manifold disconnecting of syringe and push trapped air in syringe to the waste.

If all air is not removed from the tubing, air may be forced into the patient's vascular system and cause embolism.

10. Reconnect the syringe and open the stopcocks to the syringe and repeat the procedure with loading of contrast media and saline.

Step 9 and 10 apply for set with check valve. For set without check valve continue to step 11.

11. When more bubbles appear in line or in manifold or in syringe they can be removed through tubing line by opening cap on luer lock located on end of tubing line and locking of stopcocks to drip chamber line and contrast media line.

12. Open stopcock to the patient and push mixture of fluids with syringe to the system to ensure all the manifold and lines are washed.

13. Close the stopcocks and roller clamps to close the lines.

14. Remove the protective cap from the CMD tubing luer lock.

Fig. 4

15. Connect Contrast media delivery set via luer lock to the central or peripheral catheter, *see figure 4.*

Check for fluid leakage before and during the procedure. Leaks can result in the loss of sterility, fluid and/or air embolism. If a product leaks before or during use, retighten the leaking connection or replace the product.

4. STORAGE AND DISPOSAL:

Keep dry

Use by date

Do not re-sterilise, *Improper reuse of the medical device can lead to potentially spread life-threatening infections*

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