1. GENERAL INFORMATION

Intended Use: Total Parenteral Nutrition Set is intended for the infusion, transfer, and/or administration of medications, drugs or fluids, and/or parenteral nutrition media.

Indications: is indicated when there is an inability to provide adequate nutrition via the oral / enteral route and is administered continuously until enteral feeding can be re-established.


Contraindications: Refeeding syndrome caused by rapid refeeding after a period of malnutrition, which leads to metabolic and hormonal changes characterized by electrolyte shifts (decreased phosphate, magnesium, and potassium in serum levels) that may lead to widespread cellular dysfunction.

Medical personnel responsibility for patient safety: In order to competently manage a patient, the medical personnel is required to:

- respond appropriately (adjust treatments as necessary)
- A patient on TPN must have blood level monitored closely to prevent the complications of refeeding syndrome.
- observe for complications and troubleshoot as needed

Warning notices and symbols:

The following symbols are used for important information in the Instruction for use, on the packaging:

	Caution, consult accompanying documents		Use by date
	Manufactured by		Sterilized using ethylene oxide
	CE mark of compliance		Do not use if package is damaged
	Handling mark on package – FRAGILE		For single use only
	Handling mark on package – KEEP DRY		Recyclable material
	Handling mark on package – THIS SIDE UP		Consult instruction for use
	Do not re-sterilise		Non latex
	PHT free		

2. DEVICE DESCRIPTION

TPN sets consist of components:

Tubeing, Y connector, Check valve, Needle free valve, 0,2 and 1,2 mic filter, Slide clamps, Luer lock

3. OPERATION




-  Its highly recommended to consult instruction prior to use
- 1.**  Visually check packaging integrity
- 2.**  Product is sterile, do not use if package damaged or open
- 3.**  Using antiseptic technique, open the set from the sterile packing.
- 3.** Ensure that all connections are secure, *as shown on Fig. 1.*
-  All connections should be finger tightened. Over tightening can cause cracks and leaks to occur that could result in embolism and or exposure to biohazards.
- 4.** Close the slide clamps and place the set on sterile barrier
- 5.** Remove the protective covering (pull-tab) from the TPN bag spike port and place it back on the sterile drape. Pick up the TPN set and by holding it below the sharp end of the spike tip in one hand, remove the protective covering (see Figure 2) by free hand.
-  Exposed end of drip chamber spike should not come into contact with anything, otherwise it can be contaminated.

Fig. 1

Fig. 2.

EN Instructions for use – Total Parenteral Nutrition Set

6. Insert the tip of the administration set completely into the spike port of the TPN bag using a twisting motion (see Figure 3).
7. Hang the TPN bag on an IV pole or hook (see Figure 4). In case when portable pump is used there is no need to use IV pole.
8. Squeeze the drip chamber of the administration set few times until it is half full of solution.
9. Open the slide clamp and let the solution fill the entire length of the tubing, including the filter. The filter is a safety device that traps air bubbles and any particles while the solution is infusing.
- ⚠ If all air is not removed from the tubing, air may be forced into the patient's vascular system and cause embolism.
10. When the tubing is filled, close the slide clamp.
11. Remove the protective cap from the TPN tubing luer lock.
12. Carefully attach the TPN tubing via luer lock to the central or peripheral catheter (see Figure 5).
- ⚠ Check for fluid leakage before and during the procedure. Leaks can result in the loss of sterility, fluid and/or air embolism. If a product leaks before or during use, retighten the leaking connection or replace the product.

Fig. 3.

Fig. 4.

Fig. 5.

4. STORAGE AND DISPOSAL:

-  Keep dry
-  Use by date
-  Do not re-sterilize, *Improper reuse of the medical device can lead to potentially spread life-threatening infections*

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