1. GENERAL INFORMATION

Intended Use: Connecting Lines are intended for the infusion, transfer, and/or administration of medications, drugs or fluids or parenteral nutrition media.

Indications: Broad spectrum of medical conditions where patient need to undergo infusion, transfer or administration of medications, drugs, fluids or parenteral nutrition media.

Contraindications: Not known

Medical personnel responsibility for patient safety:

In order to competently manage a patient, the medical personnel is required to:

- respond appropriately (adjust treatments as necessary)
- observe for complications and troubleshoot as needed

Warning notices and symbols: the following symbols are used for important information in the Instruction for use, on the packaging and on the product

	Caution, consult accompanying documents		Use by date
	Manufactured by		Sterilized using ethylene oxide
	CE mark of compliance		Do not use if package is damaged
	Handling mark on package – FRAGILE		For single use only
	Handling mark on package – KEEP DRY		Recyclable material
	Handling mark on package – THIS SIDE UP		Consult instruction for use
	Do not re-sterilise		Non latex
	PHT free		

2. DEVICE DESCRIPTION

The Connecting Lines consist of different lengths with various connectors and valves. The main components of the Connecting Lines are: *tubing, Luer Locks, Caps for Luer Locks*

3. OPERATION

 Its highly recommended to consult instruction prior to use.

1. Visually check packaging integrity.

 Product is sterile, do not use if package damaged or open.

2. Using antiseptic technique, open the medical device from the sterile packing

3. When connecting to accessories or other medical devices via Luer Locks ensure that all connections are secure, as shown on the Fig. 1.

Fig. 1

 All connections should be finger tightened. Over tightening can cause cracks and leaks to occur that could result in embolism and or exposure to biohazards.

 Check for fluid leakage before and during the operation. Leaks can result in the loss of sterility, fluid and/or air embolism. If a product leaks before or during use, retighten the leaking connection or replace the product.

4. STORAGE AND DISPOSAL:

 Keep dry

 Use by date

 Do not re-sterilize, *Improper reuse of the medical device can lead to potentially spread life-threatening infections*

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