1. GENERAL INFORMATION

Intended Use: *Physiology Pressure Transducer sets are intended for continuous blood pressure monitoring through incorporated pressure transducers.*

Indications:

- Induced, on-going or anticipated hypotension, or wide variations in blood pressure
- End-organ disease requiring precise pressure regulation
- The need for frequent or multiple blood gas measurements
- The need for continuous monitoring of cardiac output and stroke volume
- Situations when noninvasive methods of blood pressure monitoring are unreliable or difficult, such as with burns, trauma, or dysrhythmias

Contraindications:

- Direct blood monitoring should be avoided in extremities with suspected or pre-existing vascular insufficiency.
- Catheterization should be avoided in smaller end-arteries with inadequate collateral blood flow.

Medical personnel responsibility for patient safety:

In order to competently manage a patient the medical personnel is required to:

- respond appropriately (adjust treatments as necessary)
- understand the components of the waveform
- observe for complications and
- troubleshoot as needed

Warning notices and symbols: The following symbols are used for important information in the Instruction for use, on the packaging and on the product:

Caution, consult accompanying documents

Manufactured by

CE mark of compliance

Handling mark on package – FRAGILE

Handling mark on package – KEEP DRY

Handling mark on package – THIS SIDE UP

Do not re-sterilise

PHT free

Use by date

Sterilized using ethylene oxide

Do not use if package is damaged

For single use only

Recyclable material

Consult instruction for use

Non latex

Defibrillation-proof type CF applied part

2. DEVICE DESCRIPTION

Physiology pressure transducer sets consist of tubing sets containing disposable pressure transducer in different tubing lengths and configurations. The set can be supplied with a one channel (product codes: BP-S0xxxx, BP-SPxxxx, BP-SNxxxx), two channels (product codes: BP-D0xxxx, BP-DPxxxx) or three channel configurations (product code: BP-Txxxx) to fit customer specifications.

Main components: *3 way stopcock, Infusion/flushing system, Disposable Transducer, IV set with drip chamber, Transducer cable to signal processor, line to patient, roller clamp, Luer locks.*

3. OPERATION

Its highly recommended to consult instruction prior to use

1. Connect the adapting cable to the monitor
2. Visually check packaging integrity

Product is sterile, do not use if package damaged or open

3. Using antiseptic technique, open the kit from the sterile packing
4. Ensure that all connections are secure.

All connections should be finger tightened. Over tightening can cause cracks and leaks to occur that could result in embolism and/or exposure to biohazards.

Connect the transducer cable to the reusable cable via built hooded RJ-11 connector

- ⚠ The plug housing and the plug pins of transducer cable must not come into contact with liquid, as this can cause the corrosion of pins
- ⚠ Cable disconnection during use means wrong reading or no reading, additional treatment is needed for patient
- 6. Prime the tubing set and remove air from the IV flush solution bag per hospital policy
- ⚠ If all air is not removed from the bag, air may be forced into the patient's vascular system when the solution is exhausted
- 7. Close the roller clamp on the IV set and connect the IV set to the IV flush bag. Hang the bag approximately 60 cm above the patient. This height will provide approximately 45 mm Hg of pressure to the prime setup.
- 8. Fill the drip chamber halfway with flush solution by squeezing the drip chamber. Open the roller clamp.
- 9. To fill system:
 - a) For transducers without integral flush device fill system per hospital policy
 - b) For transducers with integral flush device, squeezing the flush device provides flow
- 10. If the kit includes IV set, open the transducer vent port by turning the stopcock handle. Deliver flush solution first through the transducer and out through the vent port, then through the set monitoring line
- ⚠ Make sure no air enters the patient through stopcock that and inadvertently left open, from accidental disconnection of the pressure setup, or from flushing residual air bubbles into the patient
- ⚠ Check for fluid leakage before and during the procedure. Leaks can result in the loss of sterility, fluid or blood loss, and/or air embolism. If a product leaks before or during use, retighten the leaking connection or replace the product

Operation with infusion syringe pump

1. Connect the *Physiology Pressure Transducer Set* to the syringe.
2. Place the syringe in the infusion syringe pump.
3. Use the infusion syringe pump in "fluid push" option.
4. If the transducer has an integral flush device, priming will be done while squeezing the flush device until solution will prime the set.
5. Replace all vented caps on side-ports of the stopcocks with non-vented caps.
6. Mount the transducer either on the patient's body per hospital procedure or to the holder.
7. Pressurize the IV flush solution bag. Flow rate will be managed by infusion syringe pump.
8. Connect pressure tubing to the catheter.
9. Flush system per hospital policy

Zeroing and Calibration

1. To ensure consistency and accuracy of the arterial blood pressure monitoring position the pressure transducer in such way that the vent port of the integrated three-way stopcock is at the level of an anatomically consistent site - phlebostatic axis.
2. Turn stop-cock on transducer so that it is off to the patient. Remove the non-vented cap and open the vent port to the atmosphere.
- ⚠ If stopcock handle is in wrong position, no reading or incorrect reading are obtained, medical procedure is not performed properly
3. Adjust the monitor to read zero mmHG.
4. Check the monitor calibration using procedure recommended by the monitor manufacturer.
5. Close the vent port to the atmosphere end replace the non-vented cap.
6. System is ready to begin motoring pressure

Blood Sampling using a syringe

1. Temporarily silence monitor alarms.
2. Unscrew the cap from luer lock on patient line and place on inside of gauze packet.
3. Connect the 3 or 5ml syringe as required
4. Turn the stopcock off to the flush bag and withdraw 3mls of blood to clear dead space in the line as required, more if dead space is bigger.
5. Turn the stopcock off to the patient and remove the syringe.
6. Connect the 10 or 20ml syringe, depending on the volume of blood required.
7. Turn the stopcock off to the flush bag and gently withdraw the required amount of blood.
8. Turn the stopcock off to the open port and remove the syringe.

9. Flush the line clear of blood using the transducer flush device.
10. Turn the stopcock off to the patient and flush blood out of the port.
11. Turn the stopcock off to the open port and check that the arterial trace returns on the monitor. Re zero the arterial line to the monitor if required.
12. Reapply or replace the cap to the port as appropriate

Testing Dynamic Responses

If there is suspicion that the values displayed on the bedside monitor are inaccurate or if the tracings are not clear one can verify the dynamic response (accuracy) of the pressure system by performing the Fast Flush test or Square Wave test.

1. Squeeze the pull or snap tab on the flush device rapidly then quickly release it
2. The monitor should show a waveform that rises suddenly and sharply, tops off, then declines sharply. As you release, one or two oscillations appear above and below the baseline after release. evaluate dynamic response

Poor dynamic response can be caused by air bubbles, clotting, excessive length of tubing, excessively compliant pressure tubing. Small bore tubing, loose connections, or leaks.

4. TROUBLESHOOTING

PROBLEM	POSSIBLE CAUSE(S)	SOLUTION/MANAGEMENT
Difficulty with zeroing - Does not reach 0 waveform - Does not reach baseline - Zeroing but not getting the number display on the monitor	Not connected properly. Line not primed appropriately through the transducer. Faulty transducer or pressure cable. Number display not set up/selected on the monitor for ABP. System not open to air	Check all equipment and connections with monitor. Ensure all roller clamps are open. Check system for air bubbles and blood clots. Recalibrate. Replace transducer, cable and/or module Replace arterial line (last resource)
Blockage/Clotting/Air embolism	Kinked line. Disconnection or loose connection. Not enough pressure on the pressure bag. Insufficient fluid in the 0.9% Sodium Chloride on the pressure bag. Line not flushed when there is the presence of blood in the line. Line not primed properly.	To prevent: Keep pressure bag inflated to ensure a pressure of 300mmHg (3-5ml/h). Ensure all connections are secure. Maintain line and blood collection port free of blood at all times. Use fast flush device to clear line to prevent clot formation. Never flush a clot or air embolism into the patient. If clot/air embolism: Attempt to aspirate blood/air to remove clot.
Falsely high or low readings Dampened wave (Over or under dampened)	Incorrect placement or transducer, either too high or too low in relation to the patient's phlebostatic axis. Uncalibrated system.	Perform square waveform check. Check position of transducer. Re-zero. Re position patient's arm. Apply back slab to maintain extension on entry point. Remove kink. Unclamp rolled clamp(s) in line Remove air bubbles/ blood clots. Check adequate amount of fluid in the flush bag. Ensure pressure bag is inflated to 300mmHg.

5. STORAGE AND DISPOSAL:

-  Keep dry
-  Use by date
-  Do not re-sterilize, *Improper reuse of the medical device can lead to potentially spread life-threatening infections*
-  Invasive blood pressure monitoring sets shall be changed every 72 hours

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