1. GENERAL INFORMATION

Intended Use: Transfer Sets Fluid Set is intended for the infusion, transfer, and/or administration of medications, drugs or fluids.


Indications: Patient medical condition in needs of infusion.

Contraindications: Not known

Medical personnel responsibility for patient safety: In order to competently manage a patient the medical personnel is required to:

- respond appropriately (adjust treatments as necessary)
- observe for complications and troubleshoot as needed

Warning notices and symbols: The following symbols are used for important information in the Instruction for use, on the packaging and on the product:

	Caution, consult accompanying documents		Use by date
	Manufactured by		Sterilized using ethylene oxide
	CE mark of compliance		Do not use if package is damaged
	Handling mark on package – FRAGILE		For single use only
	Handling mark on package – KEEP DRY		Recyclable material
	Handling mark on package – THIS SIDE UP		Consult instruction for use
	Do not re-sterilize		Non latex
	PHT free		

2. DEVICE DESCRIPTION

Transfer Set Fluids may consist of various components: *Vented Spike, Needle free valve, Y connector, Clamp, Luer Locks, Tubing, Cap for Luer Locks.*

3. OPERATION

 Its highly recommended to consult instruction prior to use

1. Visually check packaging integrity

 Product is sterile, do not use if package damaged or open

2. Using antiseptic technique, open the set from the sterile packing.

3. Ensure that all connections are secure, *as shown on Fig. 1.*

 All connections should be finger tightened. Over tightening can cause cracks and leaks to occur that could result in embolism and or exposure to biohazards.

4. Remove protective covering from spike tip and take special attention that exposed end of spike will not come into contact with anything, per Fig. 2

 Exposed end of drip chamber spike should not come into contact with anything, otherwise it can be contaminated.

5. Connect non vented bottle with medication, drugs to spike set

6. Connect syringe to needle free valve.

7. Prime the tubing line per hospital policy.

8. Evacuate the air from the tubing line.

 If all air is not removed from the tubing, air may be forced into the patient's vascular system and cause embolism

 Check for fluid leakage before and during the procedure. Leaks can result in the loss of sterility, fluid and/or air embolism. If a product leaks before or during use, retighten the leaking connection or replace the product

9. Connect tubing to the catheter via luer lock, *as shown on Fig. 3.*

10. Flush system again per hospital policy

4. STORAGE AND DISPOSAL:

 Keep dry

 Use by date

 Do not re-sterilize, *Improper reuse of the medical device can lead to potentially spread life-threatening infections*

Fig. 1

Fig. 2


Fig. 3

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