en	Instructions for use	Contrast Media Delivery Set	2
de	Gebrauchsanweisung	Kontrastmittelinjektionsset	5
fr	Mode d'emploi	Kit d'administration de produit de contraste	8
es	Instrucciones de uso	Equipo de administración de medios de contraste	11
it	Istruzioni per l'uso	Set di somministrazione di mezzi di contrasto	14
el	Οδηγίες χρήσης	Σετ χορήγησης σκιαγραφικού	17
cs	Návod k použití	Souprava pro dodávání kontrastních médií	20
sk	Návod na použitie	Súprava na podanie kontrastnej látky	23

1. GENERAL INFORMATION

Intended purpose: Contrast Media Delivery Set is intended for the infusion and/or administration of contrast media. Specifically, for injection of contrast media to patients.

Indications: It is indicated for a patient that undergoes computed tomography (contrast CT), projectional radiography, and fluoroscopy procedures, in aim to present internal organs by minimal invasive approach.

Contraindications: The risks associated to performing peripheral arteriography can be attributed to cardiac and non-cardiac complications. Specific disease states pertaining to the general medical profile of the patient (older age, renal insufficiency, uncontrolled diabetes mellitus, and morbid obesity) can increase the risk for complications. The underlying vascular status of the patient can further predispose to adverse events.

Patient target groups: general patient population.

Intended users and environment: Various hospital departments such as intensive care, emergency rooms, and diagnostic units. The use of the device should be limited to the specialist trained to perform the procedures for which this device is intended.

Medical personnel responsibility for patient safety: In order to competently manage a patient, the medical personnel are required to:

- respond appropriately (adjust treatments as necessary)
- observe for complications and troubleshoot as needed

Symbols: The following symbols are used for important information in the Instruction for use, on the packaging and on the product:

	Caution		Consult instructions for use
	Medical device		DEHP free
	Manufacturer		Does not contain, and non-presence of natural rubber latex
	Date and country of manufacture		Single sterile barrier system. Sterilized using ethylene oxide
	CE mark of compliance		Do not re-sterilize
	Use-by date		Do not use if package is damaged
	Batch code		Keep away from sunlight
	Catalogue number		Keep dry
	Unique device identifier		Handling mark on package – Fragile
	Non-pyrogenic		Handling mark on package – Recyclable material
	Do not re-use		Handling mark on package – This side up

2. DEVICE DESCRIPTION

Contrast Media delivery Set may consist of various components: Connecting Lines, Roller Clamp, Syringe, Luer Locks, Spike Set, High Pressure Manifold, Contrast Media Set, Dual Check Valve, Drip Chamber Set.

3. WARNINGS AND PRECAUTIONS

- The device is designed and intended for single use only. Do not autoclave, reshape, re-sterilize or reuse. Exposure to high heat may cause changes to the functionality of the device.

- Do not re-sterilize. Improper reuse of the medical device can lead to the potential spreading of life-threatening infections.
- Use an aseptic technique and follow hospital standard operating procedures when handling the product.
- Visually inspect the disinfected surface for possible cracks due to exposure to alcohol.
- Verify the device material is suitable for all drugs that will be used. This device contains PVC.
- This device is not intended for blood or blood products or enteral solutions.
- When using lipids, stopcocks/manifolds should be replaced every 24 hours (or per institutional protocol).
- Serious adverse events should be reported to Sky Medical and the applicable Competent Authority.

4. OPERATION

It is highly recommended to consult the instructions prior to use.

1. Visually check packaging integrity.
Product is sterile, do not use if package is damaged or open.

2. Using antiseptic technique, open the set from the sterile packing.
 3. Ensure that all connections are secure, as shown in the Fig. 1.
- All connections should be finger tightened. Over tightening can cause cracks and leaks to occur that could result in embolism and or exposure to biohazards.
4. Remove protective covering from the spike tip and take special care that the exposed end of the drip chamber spike will not come into contact with anything. Repeat the same procedure with the Drip Chamber set spike (see Figure 2).

The exposed end of the drip chamber spike should not come into contact with anything, otherwise it can be contaminated.

5. Connect the contrast line to the spike set, the saline bag to the fluid reservoir drip chamber set (see Figure 3).
6. Connect the syringe to the dual check valve through the luer lock connection.
7. Clear the contrast line and tubing line with saline by turning the stopcock on the manifold unit so both lines are opened to the syringe and locked to the patient line.
8. Evacuate the air from the contrast media set and tubing line with saline by using the syringe.
9. Evacuate the air from the syringe by closing the stop cock on the manifold, disconnecting the syringe and push the trapped air in the syringe to the waste.

If all air is not removed from the tubing, air may be forced into the patient's vascular system and cause an embolism.

10. Reconnect the syringe and open the stopcocks to the syringe and repeat the procedure with loading of the contrast media and saline.

Step 9 and 10 apply for the set with the check valve. For the set without the check valve continue to step 11.

11. When more bubbles appear in the line or in the manifold or in the syringe, they can be removed through the tubing line by opening the cap on the luer lock located on the end of the tubing line and locking the stopcocks to the drip chamber line and the contrast media line.
12. Open the stopcock to the patient and push the mixture of fluids with the syringe to the system to ensure all the manifold and lines are washed.
13. Close the stopcocks and roller clamps to close the lines.
14. Remove the protective cap from the CMD tubing luer lock.


Fig. 1

Fig. 2

Fig. 3

Fig. 4

15. Connect the Contrast media delivery set via luer lock to the central or peripheral catheter (see Figure 4).

Check for fluid leakage before and during the procedure. Leaks can result in the loss of sterility, fluid and/or air embolism. If a product leaks before or during use, retighten the leaking connection or replace the product.

5. STORAGE AND DISPOSAL

Keep dry

Use by date

After having contact with blood, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

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