

RondelO™ End-line adapter

EN	Instructions for Use	2	RO	Instrucțiuni de utilizare	62
DE	Gebrauchsanweisung	6	SL	Navodila za uporabo	66
ES	Instrucciones de uso	10	SV	Bruksanvisning	70
FR	Mode d'emploi	14	CS	Návod k použití	74
IT	Istruzioni per l'uso	18	HU	Használati útmutató	78
NL	Gebruiksaanwijzing	22			
PT	Instruções de utilização	26			
FI	Käyttöohjeet	30			
SK	Návod na použitie	34			
IL	שומי של תוארה	38			
EL	Οδηγίες χρήσης	42			
DA	Bruksanvisning	46			
HR	Upute za uporabu	50			
NO	Bruksanvisning	54			
PL	Instrukcja obsługi	58			

1. General Information

Intended purpose: The RondelO End-line adapter is intended for the infusion/administration of various commonly used medical fluids.

Indications: Patient medical condition in need of infusion.

Contraindications: Once an intravenous medication is delivered, it cannot be retrieved. When giving IV medications, there is very little opportunity to stop an injection if an adverse reaction or error occurs. IV medications, if given too quickly or incorrectly, can cause significant harm.

Patient target groups: General patient population

Intended users and environment: Various hospital departments such as intensive care, emergency rooms and diagnostic units. The use of the device should be limited to the specialist trained to perform the procedures for which this device is intended.

Medical personnel responsibility for patient safety: In order to competently manage a patient, the medical personnel are required to:

- respond appropriately (adjust treatments as necessary)
- observe for complications and troubleshoot as needed

Symbols: The following symbols are used for important information in the Instructions for Use, on the packaging and on the product.

	Caution		Consult Instructions for Use
	Medical device		DEHP free
	Manufacturer		Does not contain, and there is no presence of, natural rubber latex
	Date and country of manufacture		Do not re-use
	CE mark of compliance		Single sterile barrier system. Sterilized using ethylene oxide
	Use-by date		Do not re-sterilize
	Batch code		Do not use if package is damaged
	Catalogue number		Keep away from sunlight
	Unique device identifier		Keep dry
	Non-pyrogenic		Handling mark on package – Fragile
	Drops per millilitre		Handling mark on package - Recyclable material
	Liquid filter with 15-micron pore size		Handling mark on package – This side up
	One-way valve		

Abbreviations

DC - Drip Chamber

FLL - Female Luer Lock

IV-set - Intravenous infusion set

NFV - Needle Free Valve

MLL - Male Luer Lock

2. Warnings and Precautions



- The device is designed and intended for single use only. Do not autoclave, reshape, re-sterilize or reuse. Exposure to high heat may cause changes to the functionality of the device.
- Do not re-sterilize. Improper reuse of the medical device can lead to the potential spreading of life-threatening infections.
- Sets are to be changed every 72 hours.
- Make sure the primary infusion is suitable as a flushing fluid for all drugs that will be used.
- Verify the device material is suitable for all drugs that will be used. This device contains PVC.
- Use aseptic technique and follow hospital standard operating procedures when handling the product.
- Respond appropriately and adjust treatments, as necessary.
- Observe for complications and troubleshoot as needed.
- Visually inspect the disinfected surface for possible cracks due to exposure to alcohol.
- Serious adverse events should be reported to Sky Medical a.s. and the applicable Competent Authority.
- This device is not intended for blood or blood products or enteral solutions.

RondelO positions

Turn the handle to select the intended position. Each position is verified by a haptic click.

- **Drug positions** – Selectable Drug positions correspond to Luer lock connectors of the RondelO device. When one Drug position is selected, the RondelO is open at the corresponding Luer lock connector, and is closed at all other Luer lock connectors.
- **Flushing positions** – Selectable Flushing positions are located in between the Drug positions to allow primary infusion to flush solely and rinse the RondelO internally.

Note: Primary infusion is available in all RondelO positions. Concurrent flushing from Primary infusion is available by default in all Drug positions. Availability of Primary infusion is controlled by the Clamp above the RondelO.



It is highly recommended to consult the instructions prior to use.

3. RondelO End-line adapter

DEVICE DESCRIPTION: RondelO End-line adapter may consist of various components:

(1) Female Luer Lock (FLL) with one-way valve, (2) Protective cap, (3) Slide clamp, (4) RondelO with 2, 4, or 6 Drug positions at corresponding Luer Lock connectors, with either a (5) Protective cap, or (6) Needle free valve (NFV), (7) RondelO Handle, (8) Pinch clamp, (9) Male Luer Lock (MLL), (10) Prime stop cap, Tubing, Transparent or UV-protected components. *Note: Drawings are not to scale.*

4. Operation

RondelO End-line adapter - priming

1. Visually check packaging integrity.
-  2. Product is sterile, do not use if package damaged or open.
3. Using aseptic technique, open the sterile packaging to access the RondelO End-line adapter.
4. Verify the RondelO Handle (7) is in a Flushing position between two Drug positions.
5. Remove the Protective cap (2) from the Female Luer Lock with a one-way valve (1) on the RondelO End line adapter.
-  The exposed end of the Female Luer Lock should not come into contact with anything, otherwise it can be contaminated.
6. Connect a Primary infusion to the Female Luer Lock with a one-way valve (1).
-  All connections should be finger tightened. Over tightening can cause cracks and leaks to occur that could result in an air embolism and or exposure to biohazards.
7. Verify or open the slide clamp (3) and the pinch clamp (8) and prime the RondelO End-line adapter.
8. Close the Pinch clamp (8) of the RondelO end-line adapter.
9. Remove the Prime stop cap (10) from the Male Luer Lock (9) of the RondelO End line adapter.
-  If all air is not removed from the tubing, air may be forced into the patient's vascular system and cause an air embolism.
10. Connect the end-line adapter to the patient's IV access port.
-  Check for fluid leakage before and during the procedure. Leaks can result in the loss of sterility, fluid and/or air embolism. If a product leaks before or during use, retighten the leaking connection or replace the product

RondelO End-line adapter – connect a secondary IV-set or a syringe to the Drug position to be used of the RondelO

1. Verify that the RondelO Handle (7) is not aligned with the Drug position.
2. If applicable, remove the protective cap (5) from the Luer Lock connector at the Drug position.
3. Connect the secondary IV-set or syringe to the Luer lock connector at the Drug position.

RondelO End-line adapter – administration

1. Close or open the Slide clamp (3) on the RondelO End-line adapter depending on whether the drug at the selected Drug position should be administered with concomitant flushing or not.
2. Turn the RondelO Handle (7) to the Drug position.
3. Open the Pinch clamp (8) towards the patient.
4. Administer the drug.

Flushing between and after Drug administrations

1. Turn the RondelO Handle (7) to an adjacent Flushing position.
2. Open the Slide clamp (3) of the RondelO End-line adapter and flush.

Disconnect from the Luer Lock connector at a Drug position of the RondelO

1. Verify that the RondelO Handle (7) is not aligned with the Drug position to be disconnected.
2. If applicable, verify the clamp on the secondary IV-set is closed, before disconnection.
3. Disconnect the secondary IV-set or syringe from RondelO End-line adapter.
4. If applicable, cap the Luer lock connector at the drug position of RondelO.

Disconnect RondelO End-line adapter from patient

1. Close the Pinch clamp (8) of the RondelO End-line adapter.
2. Disconnect the RondelO End-line adapter from the patient's IV access port.

5. Storage and Disposal

Keep dry

Handling mark on package – Fragile

Use by date

Keep away from sunlight

After having contact with the blood, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

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