
RondelO™ Primary IV-set

EN	Instructions for Use	2	RO	Instrucțiuni de utilizare	77
DE	Gebrauchsanweisung	7	SL	Navodila za uporabo	82
ES	Instrucciones de uso	12	SV	Bruksanvisning	87
FR	Mode d'emploi	17	CS	Návod k použití	92
IT	Istruzioni per l'uso	22	HU	Használati útmutató	97
NL	Gebruiksaanwijzing	27			
PT	Instruções de utilização	32			
FI	Käyttöohjeet	37			
SK	Návod na použitie	42			
IL	הוראות שימוש	47			
EL	Οδηγίες χρήσης	52			
DA	Bruksanvisning	57			
HR	Upute za uporabu	62			
NO	Bruksanvisning	67			
PL	Instrukcja obsługi	72			

1. General Information

Intended purpose: The RondelO Primary IV-set is intended for an infusion/administration of various commonly used medical fluids.

Indications: Patient medical condition in need of infusion.

Contraindications: Once an intravenous medication is delivered, it cannot be retrieved. When giving IV medications, there is very little opportunity to stop an injection if an adverse reaction or error occurs. IV medications, if given too quickly or incorrectly, can cause significant harm.

Patient target groups: General patient population

Intended users and environment: Various hospital departments such as intensive care, emergency rooms and diagnostic units. The use of the device should be limited to the specialist trained to perform the procedures for which this device is intended.

Medical personnel responsibility for patient safety: In order to competently manage a patient, the medical personnel are required to:

- respond appropriately (adjust treatments as necessary)
- observe for complications and troubleshoot as needed

Symbols: The following symbols are used for important information in the Instructions for Use, on the packaging and on the product.

	Caution		Consult Instructions for Use
	Medical device		DEHP free
	Manufacturer		Does not contain, and there is no presence of, natural rubber latex
	Date and country of manufacture		Do not re-use
	CE mark of compliance		Single sterile barrier system. Sterilized using ethylene oxide
	Use-by date		Do not re-sterilize
	Batch code		Do not use if package is damaged
	Catalogue number		Keep away from sunlight
	Unique device identifier		Keep dry
	Non-pyrogenic		Handling mark on package – Fragile
	Drops per millilitre		Handling mark on package - Recyclable material
	Liquid filter with 15-micron pore size		Handling mark on package – This side up
	One-way valve		

Abbreviations

DC - Drip Chamber

FLL - Female Luer Lock

IV-set - Intravenous infusion set

NFV - Needle Free Valve

MLL - Male Luer Lock

2. Warnings and Precautions


- The device is designed and intended for single use only. Do not autoclave, reshape, re-sterilize or reuse. Exposure to high heat may cause changes to the functionality of the device.
- Do not re-sterilize. Improper reuse of the medical device can lead to the potential spreading of life-threatening infections.
- Sets are to be changed every 72 hours.
- Make sure the primary infusion is suitable as a flushing fluid for all drugs that will be used.
- Verify the device material is suitable for all drugs that will be used. This device contains PVC.
- Use aseptic technique and follow hospital standard operating procedures when handling the product.
- Respond appropriately and adjust treatments, as necessary.
- Observe for complications and troubleshoot as needed.
- Visually inspect the disinfected surface for possible cracks due to exposure to alcohol.
- Serious adverse events should be reported to Sky Medical a.s. and the applicable Competent Authority.
- This device is not intended for blood or blood products or enteral solutions.

RondelO positions

Turn the handle to select the intended position. Each position is verified by a haptic click.

- **Drug positions** – Selectable Drug positions correspond to Luer lock connectors of the RondelO device. When one Drug position is selected, the RondelO is open at the corresponding Luer lock connector, and is closed at all other Luer lock connectors.
- **Flushing positions** – Selectable Flushing positions are located in between the Drug positions to allow primary infusion to solely flush and rinse the RondelO internally.


Note: Primary infusion is available in all RondelO positions. Concurrent flushing from Primary infusion is available by default in all Drug positions. Availability of Primary infusion is controlled by the Clamp above the RondelO or head height differential.

 It is highly recommended to consult the instructions prior to use.

3. RondelO Primary IV-set (DC in line)

DEVICE DESCRIPTION: *The RondelO Primary IV-set (DC in line)* may consist of various components: (1) Vented Spike with protective covering, (2) Pinch clamp, (3) One-way valve, (4) RondelO with 2, 4 or 6, Drug positions at a corresponding Luer lock connector with either a (5) Protective cap, or (6) Needle free valve (NFV), (7) RondelO Handle, (8) Drip Chamber, (9) Roller clamp, (10) Y-connector, (11) Slide clamp, (12) Male Luer Lock (MLL), (13) Prime stop Cap, Tubing, Transparent or UV protected components. *Components 10 and 11 are not present in some of the IV-set configurations. Note: Drawings are not to scale.*

4. Operation

RondelO Primary IV-set (DC in line) – priming

1. Visually check packaging integrity.
-  Product is sterile, do not use if package damaged or open.
2. Using aseptic technique, open the sterile packaging to access the RondelO Primary IV-set.
3. When applicable close the vented cap on the spike (1).
4. Verify the RondelO handle (7) is in a Flushing position.
5. Close the Roller clamp (9) of the RondelO Primary IV-set.
6. Remove the protective covering from spike tip.
-  The exposed spike tip should not come into contact with anything, otherwise it can become contaminated.
7. Spike the Infusion bag/container.
8. Hold the Drip chamber (8) tilted upside down and position it at the same level as the RondelO, open the Roller clamp (9) of the IV-set and fill the Drip chamber as in (Fig. 1.)


Fig 1. Priming of the drip chamber

9. Continue to prime the rest of the Primary IV-set.
-  If all air is not removed from the tubing, air may be forced into the patient's vascular system and cause an air embolism.
10. Close the Roller clamp (8).
11. Remove the Prime stop cap (13) of the RondelO Primary IV-set and connect it (12) to the patient's IV-access port.
-  Check for fluid leakage before and during the procedure. Leaks can result in the loss of sterility, fluid and/or air embolism. If a product leaks before or during use, retighten the leaking connection or replace the product
-  All connections should be finger tightened. Over tightening can cause cracks and leaks to occur that could result in an air embolism and or exposure to biohazards.

RondelO Primary IV- set (DC in line) – connect a secondary IV-set to the Drug position to be used on the RondelO

1. Verify the clamp on the secondary IV-set is closed before connecting.
2. Verify the RondelO handle (7) is not aligned with the Drug position.
3. If applicable, remove the protective cap (5) from the Luer lock connector at the Drug position.
4. Remove the cap from the secondary IV-set and connect it to the Luer lock connector at the Drug position.
5. Repeat the procedure for each successive secondary IV-set.

RondelO Primary IV- set (DC in line) – administer secondary infusions

1. Close the Pinch clamp (2) on the Primary IV-set above the RondelO.
2. Turn the RondelO handle (7) to the Drug position and open the clamp on the secondary IV-set.
3. Adjust the Roller clamp (9) on the RondelO Primary IV-set to start and regulate infusion to patient.
4. Close the clamp on the secondary IV-set after the drug administration is completed.

Flushing between and after drug administrations

1. Turn the RondelO handle (7) to an adjacent Flushing position.
2. Open the Pinch clamp (2) above the RondelO and flush.

Disconnect from the RondelO


1. Verify that the RondelO handle (7) is not aligned with the Drug position to be disconnected.
2. Verify the clamp on the secondary IV-set is closed, before disconnection.
3. Disconnect the secondary IV-set from the RondelO.
4. If applicable, cap the Luer lock connector at the drug position of RondelO.

Disconnect the RondelO Primary IV-set from the patient

1. Close the Roller clamp (9) on the Primary IV-set.
2. Disconnect the IV-set from the patient's IV access port.

5. RondelO Primary IV-set (DC on top)

DEVICE DESCRIPTION: RondelO Primary IV-set (DC on top) may consist of various components: (1) Vented spike with a protective covering integrated with a Drip chamber, (2) One-way valve, (3) RondelO with 2, 4, or 6, Drug position at corresponding Luer Lock connectors, with either a (4) Protective cap, or (5) Needle free valve (NFV), (6) RondelO Handle, (7) Roller clamp, (8) Y-connector, (9) Slide clamp, (10) Male Luer Lock (MLL), (11) Prime stop Cap, Tubing. Transparent or UV protected components. *Components 8 and 9 are not present in some of the IV-set configurations. Note: Drawings are not to scale.*

6. Operation

RondelO Primary IV-set (DC on top) – priming

1. Open the sterile packaging to access the RondelO Primary IV-set.
2.  Product is sterile, do not use if package damaged or open.
3. When applicable close the vented cap on the spike (1).
4. Verify the RondelO handle (6) is in a Flushing position.
5. Close the Roller clamp (7) of the RondelO Primary IV-set.
6.  Remove the protective covering from spike tip.
The exposed spike tip should not come into contact with anything, otherwise it can become contaminated.
7. Spike the Infusion bag/container.
8. Prime the drip chamber (1) with the squeeze and release technique.
9. Open the Roller clamp (7) to prime the RondelO and the rest of the IV-set.
10.  If all air is not removed from the tubing, air may be forced into the patient's vascular system and cause an air embolism.
11. Close the Roller clamp (7).

10. Remove the prime stop cap (11) of the primary IV-set and connect it (10) to the patient's IV-access port.

Check for fluid leakage before and during the procedure. Leaks can result in the loss of sterility, fluid and/or air embolism. If a product leaks before or during use, retighten the leaking connection or replace the product.

All connections should be finger tightened. Over tightening can cause cracks and leaks to occur that could result in an air embolism and or exposure to biohazards.

RondelO Primary IV- set (DC on top) – connecting secondary IV-set to a Drug position to be used with the RondelO

1. Verify the clamp on the secondary IV-set is closed before connecting.
2. Verify the RondelO handle (6) is not aligned with the Drug position.
3. If applicable, remove the protective cap (4) from the Luer Lock connector at the Drug position.
4. Remove the cap from the secondary IV-set and connect it to the Luer Lock connector at the Drug position.
5. Repeat the procedure for each successive secondary IV-set.

RondelO Primary IV- set (DC on top) – administer secondary infusions with automatic flushing take over

1. Utilise the head height differential practise to regulate flow from primary or secondary infusion.
2. Turn the RondelO handle (6) to the Drug position and open the clamp fully on the secondary IV-set.
3. Regulate the flow rate of the infusion to the patient by adjusting the roller clamp (7) of the RondelO Primary IV-set.

Flushing fluid from the primary bag will automatically take over after the Secondary IV-infusion is completed.

Flushing between and after drug administrations

1. Turn the RondelO handle (6) to an adjacent Flushing position.

Disconnect from the RondelO

1. Verify that the RondelO handle (6) is not aligned with the Drug position to be disconnected.
2. Verify the clamp on the secondary IV-set is closed, before disconnection.
3. Disconnect the secondary IV-set from the RondelO Primary IV-set.
4. If applicable, cap the Luer lock connector at the drug position of RondelO.

Disconnect the RondelO Primary IV-set from the patient

1. Close the Roller clamp (7) on the primary IV-set.
2. Disconnect the IV-set from the patient's IV access port.

7. Storage and Disposal

Keep dry

Handling mark on package – Fragile

Use by date

Keep away from sunlight

After having contact with the blood, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

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